





Acylation of a Support-Bound Amine with an Aryl Isocyanate

Acylation of primary or secondary amines with isocyanates gives rise to substituted ureas. In the following example, a support-bound aniline¹ on SynPhase™ PS Lanterns² is acylated with an aryl isocyanate, to form a diaryl urea.



DMSO: dimethylsulfoxide
DMF: dimethylformamide
TFA: trifluoroacetic acid
DCM: dichloromethane



Urea Formation

Each D-Series Lantern derivatized with an aniline¹ (initial specified loading: 34 μmol) is treated with 0.5 mL of a 0.4 M solution of 4-trifluoromethylphenylisocyanate (220 μmol , 6.5 mole equivalents) in anhydrous DCM at

35°C for 18 h. The reaction is allowed to cool to room temperature and the solution decanted. The Lanterns are washed with DCM (1 \times 3 min), DMF (3 \times 3 min) and DCM (3 \times 3 min) then air dried.



Cleavage

Individual Lanterns are placed in polypropylene tubes and treated with 20% TFA/DCM (0.6–0.8 mL) for 1 h. The Lanterns are removed and the cleaved products are

concentrated using a centrifugal evaporator. Samples are dissolved in 90% $\text{CH}_3\text{CN}/\text{H}_2\text{O}$ for HPLC and ES-MS analysis.


Reverse-phase HPLC traces of the crude urea

Top trace: detection at 214nm
Bottom trace: detection at 254nm



Electrospray MS trace of the crude urea

Molecular Formula: $C_{15}H_{11}ClF_3N_3O_2$
Monoisotopic Mol. Weight: 357.7amu
 $[M+H]^+ = 358.8\text{amu}$
 $[2M+H]^+ = 716.4\text{amu}$



400MHz ^1H NMR spectrum of the crude urea (D_6 -DMSO)

(\otimes) DMSO
(\times) Residual DMF

References and Notes

- 1 See SynPhase Chemistry Note SCN 001-3.
- 2 The chemistry described here was performed using SynPhase PS Lanterns but is readily adaptable to SynPhase PA Lanterns.



MIMOTOPES

International
Tel: +61 3 9565 1111
Fax: +61 3 9565 1199
mimotopes@mimotopes.com

France
Tel: +33 1 5858 0002
Fax: +33 1 5858 0006
europe@mimotopes.com

United Kingdom
Tel: +44 151 648 3343
Fax: +44 151 648 3328
uk@mimotopes.com

USA West
Tel: +1 858 558 5800
Fax: +1 858 558 5810
Tel: 800 644 1866
Fax: 800 655 1866
uswest@mimotopes.com

USA East
Tel: +1 919 873 1123
Fax: +1 919 873 1127
Tel: 800 633 8161
Fax: 800 424 3970
useast@mimotopes.com