





## Preparation of a Bicyclic Pyrrolidine via a Three Component 1,3-Dipolar Cycloaddition

**Pyrrolidines can be generated** by trapping an azomethine ylide with a dipolarophile. In the following example, a three component 1,3-dipolar cycloaddition involving a solid-bound benzaldehyde,<sup>1</sup> an amino acid methyl ester and an *N*-substituted maleimide are used to prepare bicyclic pyrrolidines<sup>2</sup> on SynPhase™ PS HMP derivatized Lanterns.<sup>3</sup>



DMF: dimethylformamide  
DCM: dichloromethane  
TFA: trifluoroacetic acid  
TEA: triethylamine  
HMP: hydroxymethylphenoxy



### Three-Component 1,3-Dipolar Cycloaddition

Each **D-Series Lantern** derivatized with 4-hydroxybenzaldehyde<sup>1</sup> (initial specified loading: 36 $\mu$ mol) is treated with 0.5mL of a solution of L-leucine methyl ester hydrochloride (0.2 M, 100 $\mu$ mol, 2.8 mole equivalents), *N*-phenyl maleimide (0.2M, 100 $\mu$ mol, 2.8 mole equivalents), TEA (0.2M,


100 $\mu$ mol, 2.8 mole equivalents) and acetic acid (0.5M, 250 $\mu$ mol, 7.0 mole equivalents) in DMF at 100°C for 24h. The reaction is allowed to cool to room temperature, the reagent solution is decanted and the Lanterns washed in turn with DMF (3x3min) and DCM (3x3min) then air dried.

## Cleavage

**Individual Lanterns are placed** in polypropylene tubes and treated with 20% TFA/DCM (0.6-0.8mL) for 1h. Cleavage solutions are concentrated using a centrifugal evaporator. The crude yield is ca. 70%, based on the initial specified loading of the


Lanterns. LC-MS analysis indicates three LC peaks with the required  $[M+H]^+$  peak in the mass spectra. The purity of the crude diastereomeric mixture is 95%. Samples are dissolved in 95%  $CH_3CN/H_2O$  for LC-MS analysis.

## Analytical Data




### Electrospray Mass spectra of LC peaks at $R_t = 7.6, 8.1$ and $8.3$ min

Molecular Formula:  $C_{24}H_{26}N_2O_5$   
Monoisotopic Mol. Weight: 422.5amu  
 $[M+H]^+$  peak at 423.4amu in each spectrum



### 400MHz $^1H$ NMR spectrum of crude pyrrolidine ( $CDCl_3$ )



### References

- 1 See SynPhase Chemistry Note SCN 005-3.
- 2 Hamper, B.C., Dukeshner, D.R. and South, M.S., *Tetrahedron Lett.*, 1996, **37**, 3671-3674.
- 3 The chemistry described here was performed using SynPhase PS Lanterns but is readily adaptable to SynPhase PA Lanterns.



**International**  
Tel: + 61 3 9565 1111  
Fax: + 61 3 9565 1199  
mimotopes@mimotopes.com

**France**  
Tel: + 33 1 5858 0002  
Fax: + 33 1 5858 0006  
europe@mimotopes.com

**United Kingdom**  
Tel: +44 151 648 3343  
Fax: +44 151 648 3328  
uk@mimotopes.com

**USA West**  
Tel: + 1 858 558 5800  
Fax: + 1 858 558 5810  
Tel: 800 644 1866  
Fax: 800 655 1866  
uswest@mimotopes.com

**USA East**  
Tel: + 1 919 873 1123  
Fax: + 1 919 873 1127  
Tel: 800 633 8161  
Fax: 800 424 3970  
useast@mimotopes.com