





Preparation of Secondary Carboxamides on BAL Linker

A range of secondary carboxamides may be generated using the backbone amide forming Linker in a two step procedure.¹ In the following example, 2-aminophenol is reductively aminated to the linker on SynPhase™ PS Lanterns.² The resultant support-bound secondary amine is then acylated with a carboxylic acid such as 2-fluoro-5-nitrobenzoic acid using tetrafluoromethylformamidinium hexafluorophosphate.



BAL: backbone amide linker
DMF: dimethylformamide
DCM: dichloromethane
rt: room temperature
DIEA: diisopropylethylamine
DCM: dichloromethane
TFA: trifluoroacetic acid



Reductive Amination

Each D-series Lantern (initial specified loading: $36\mu\text{mol}$) is treated with 0.5mL of a solution of 2-aminophenol (0.5M, $250\mu\text{mol}$, 7 mole equivalents) and sodium cyanoborohydride (0.05M, $25\mu\text{mol}$, 0.7 mole

equivalents) in 1% acetic acid/DMF at 60°C for 17h. After cooling to rt, the reagent solution is decanted and the Lanterns washed with DMF ($3\times 3\text{min}$) and DCM ($3\times 3\text{min}$).

Acylation

Each D-Series Lantern is treated with 0.5mL of a solution of 2-fluoro-5-nitrobenzoic acid (0.1M, 50 μ mol, 1.4 mole equivalents), TFFH (0.1M, 50 μ mol, 1.4 mole equivalents) and DIEA (0.2M,

100 μ mol, 2.8 mole equivalents) in 1% DMF/DCM at rt for 2h. The reagent solution is decanted, the Lanterns are washed with DMF (3 \times 3min) and DCM (3 \times 3min) then air dried.

Cleavage


Individual Lanterns are placed in polypropylene tubes and treated with 20% TFA/DCM (0.6-0.8mL) for 1h. The Lanterns are removed and the cleaved products are concentrated using a centrifugal evaporator. The residue is freeze-

dried from 90% CH₃CN/H₂O. The yield of carboxamide product is 54%, based on the initial loading. Samples are dissolved in 90% CH₃CN/H₂O for HPLC and ES-MS analysis.

Analytical Data


Reverse phase HPLC trace of the crude product

Detection at 214nm



Electrospray MS trace of LC peak at R_t=7.24min

Molecular Formula: C₁₃H₉FN₂O₄
Monoisotopic Mol. Weight: 276amu
[M+H]⁺ peak at 277.1amu
[2M+H]⁺ peak at 553.4amu



References and Notes

1 Ouyang, X., Tamayo, N. and Kiselyov, A.S., *Tetrahedron*, 1999, **55**, 2827-2834.

2 The chemistry described here was performed using SynPhase PS Lanterns but is readily adaptable to SynPhase PA Lanterns.



International
Tel: + 61 3 9565 1111
Fax: + 61 3 9565 1199
mimotopes@mimotopes.com

France
Tel: + 33 1 5858 0002
Fax: + 33 1 5858 0006
europe@mimotopes.com

United Kingdom
Tel: +44 151 648 3343
Fax: +44 151 648 3328
uk@mimotopes.com

USA West
Tel: + 1 858 558 5800
Fax: + 1 858 558 5810
Tel: 800 644 1866
Fax: 800 655 1866
uswest@mimotopes.com

USA East
Tel: + 1 919 873 1123
Fax: + 1 919 873 1127
Tel: 800 633 8161
Fax: 800 424 3970
useast@mimotopes.com