





S_NAr Reaction of an Activated Aryl Fluoride with a Primary amine

A range of aniline derivatives may be prepared by treating a solid-bound aryl fluoride with primary amines, via S_NAr displacement. In the following example, an *o*-nitroaryl fluoride on SynPhase™ hydroxymethylphenoxy Lanterns is efficiently converted to an *N*-substituted *o*-nitroaniline. The performance of this type of reaction is enhanced when an electron withdrawing group is also present.



DMF: dimethylformamide
DCM: dichloromethane
TFA: trifluoroacetic acid
DIEA: diisopropylethylamine
DMAP: 4-(dimethylamino)pyridine

DIC: diisopropylcarbodiimide
rt: room temperature



Attachment of 4-Fluoro-3-Nitrobenzoic Acid

Each D-Series HMP Lantern (initial specified loading: 35μmol) is treated with a 0.5mL solution of 4-fluoro-3-nitrobenzoic acid (0.56M, 280μmol, 8 mole equivalents), DIC (0.28M, 140μmol, 4 mole equivalents) and

DMAP (0.02M, 10μmol, 0.3 mole equivalents) in DMF. The reaction mixture is stood at rt for 2h. The reagent solution is decanted and the Lanterns washed with DMF (3×3min), DCM (3×3min) and air dried.

S_NAr Reaction

Each of the D-Series Lanterns prepared above is treated with a 0.5mL solution of *n*-propylamine (0.35M, 175μmol, 5 mole equivalents) in 5% DIEA/DMF. The reaction

mixture is stood at rt for 6h.¹ The reagent solution is then decanted and the Lanterns washed with DMF (3×3min), DCM (3×3min) and air dried.

Cleavage


Individual Lanterns are placed in polypropylene tubes and treated with 20% TFA/DCM (0.6-0.8mL) for 1h. The Lanterns are removed and the cleaved products are concentrated using a centrifugal

evaporator, then the residues are freeze-dried from 90% CH₃CN/H₂O. The yield of product was essentially quantitative. Samples are redissolved in the same solvent for HPLC and ES-MS analysis.

Analytical Data


Reverse phase HPLC trace of the crude S_NAr product

Detection at 214nm



400MHz ¹H NMR spectrum of the crude S_NAr product (D₄-methanol)

(€) CH₃OH
(X) DMF
(⊗) H₂O



Note

1 For primary anilines, overnight heating at 60°C may be required for complete conversion.



MIMOTOPES

International
Tel: + 61 3 9565 1111
Fax: + 61 3 9565 1199
mimotopes@mimotopes.com

France
Tel: + 33 1 5858 0002
Fax: + 33 1 5858 0006
europe@mimotopes.com

United Kingdom
Tel: +44 151 648 3343
Fax: +44 151 648 3328
uk@mimotopes.com

USA West
Tel: + 1 858 558 5800
Fax: + 1 858 558 5810
Tel: 800 644 1866
Fax: 800 655 1866
uswest@mimotopes.com

USA East
Tel: + 1 919 873 1123
Fax: + 1 919 873 1127
Tel: 800 633 8161
Fax: 800 424 3970
useast@mimotopes.com



<http://www.synphase.com> <http://www.mimotopes.com>