




One-Pot Synthesis of a Benzimidazole from a Solid-Bound *o*-Nitroaniline

Benzimidazoles may be generated efficiently by heating a solid-bound *o*-nitroaniline and an aldehyde in the presence of tin(II) chloride dihydrate.¹ In the following example, an *o*-nitroaniline prepared² on hydroxymethylphenoxy derivatized SynPhase™ Lanterns is converted to a benzimidazole in a one-pot procedure.



Synthesis of Benzimidazole

Each *o*-nitroaniline derivatised-HMP D-Series Lantern is treated with a 0.5mL solution of 3-(trifluoromethyl)benzaldehyde (0.15M, 75 μ mol, 2.1 mole equivalents) and tin(II) chloride dihydrate (0.75M, 375 μ mol, 10.7

mole equivalents) and heated at 60°C for 3h.³ The reagent solution is then decanted and the Lanterns are washed with DMF (3 \times 3min), 80% THF/H₂O⁴ (60°C, 2 \times 1h), MeOH (2 \times 3min) and DCM (3 \times 3min).

Cleavage



Individual Lanterns are placed in polypropylene tubes and treated with 20% TFA/DCM (0.6-0.8mL) for 1h. The Lanterns are removed and the cleaved products are concentrated using a centrifugal evaporator. The residue is freeze-

dried from 90% CH₃CN/H₂O. The benzimidazole product was obtained in ca. 85% yield, based on the initial Lantern loading. Samples are dissolved in 90% CH₃CN/H₂O for HPLC and ES-MS analysis.

Analytical Data

Reverse phase HPLC trace of the crude benzimidazole

Detection at 214nm




Electrospray MS trace of LC peak at R_t = 8.91 min

Molecular Formula: C₁₈H₁₅F₃N₂O₂
Monoisotopic Mol. Weight: 348.1amu
[M+H]⁺ peak at 349.1amu

400MHz ¹H NMR spectrum of the crude benzimidazole (CDCl₃)

(x) Residual DMF



References & Notes

- 1 Wu, Z., Rea, P. and Wickham, G., *Tetrahedron Lett.*, 2000, **41**, 9871-9874.
- 2 See Mimotopes SynPhase Chemistry Note SCN 014-2.
- 3 If the solid supported species is *N*-diphenyl substituted, then a longer reaction time may be required.
- 4 THF/H₂O washing assists in the removal of tin residues.



International
Tel: + 61 3 9565 1111
Fax: + 61 3 9565 1199
mimotopes@mimotopes.com

France
Tel: + 33 1 5858 0002
Fax: + 33 1 5858 0006
europe@mimotopes.com

United Kingdom
Tel: +44 151 648 3343
Fax: +44 151 648 3328
uk@mimotopes.com

USA West
Tel: + 1 858 558 5800
Fax: + 1 858 558 5810
Tel: 800 644 1866
Fax: 800 655 1866
uswest@mimotopes.com

USA East
Tel: + 1 919 873 1123
Fax: + 1 919 873 1127
Tel: 800 633 8161
Fax: 800 424 3970
useast@mimotopes.com