




Preparation of Tertiary Amines via Acrylate-based REM Linker on SynPhase™ PS Hydroxymethyl Lanterns

Both sulfonylate-based¹ and acrylate-based² REM linkers have been used in the solid phase synthesis of tertiary amines from secondary amines. The synthesis involves (1) coupling the starting secondary amine to the REM linker (Michael addition), (2) modifying the attached amine (3) quaternization of the resultant tertiary amine (activation) and (4) cleavage of the product with a base such as diisopropylethylamine (Hofmann elimination). A shortcoming of this synthesis is that the final tertiary amine product is contaminated with the acid salt of the base. To remove the contaminant, chromatography on silica gel or ion-exchange resin is usually required.

In the following example, a tertiary amine is prepared via the acrylate REM linker on SynPhase™ PS hydroxymethyl Lanterns using the above mentioned chemistries.

The Hofmann elimination cleavage, however, is accomplished using both conventional base² (Method A) and PS aminomethyl Lanterns as a solid phase reagent^{3,4} (Method B). Using PS aminomethyl Lanterns as a base shows great advantage since no further purification is required to obtain the pure product.




DMF:	dimethylformamide	REM:	class of acrylate based linkers
DCM:	dichloromethane	OH :	hydroxymethyl Lantern
DIEA:	diisopropylethylamine	NH ₂ :	aminomethyl Lantern
rt:	room temperature		
TFA:	trifluoroacetic acid		
NMP:	N-methyl pyrrolidone		

Preparation of REM Lanterns from hydroxymethyl Lanterns

Each PS hydroxymethyl D-Series Lantern (initial specified loading, $36\mu\text{mol}$) is treated with 0.5 mL of a solution of DIEA ($50\mu\text{L}$, $287\mu\text{mol}$, 8 mole equivalents) in DCM. To the mixture is added acryloyl chloride ($26\mu\text{L}$, $287\mu\text{mol}$, 8 mole


equivalents). The reaction is allowed to stand at rt for 4 hours. The reagent solution is decanted and the Lanterns washed with DMF ($3\times 3\text{min}$) and DCM ($3\times 3\text{min}$).



Michael Addition

Each D-Series acrylate REM Lantern is treated with 0.5 mL of a solution of piperazine ($50\mu\text{L}$, $287\mu\text{mol}$, 8 mole equivalents) in DMF at rt for 16


hours. The reagent solution is decanted and the Lanterns washed with DMF ($3\times 3\text{min}$) and DCM ($3\times 3\text{min}$).



Alkylation

Each D-Series Lantern is treated with 0.5 mL of a solution of DIEA ($62\mu\text{L}$, $360\mu\text{mol}$, 10 mole equivalents), and bromodiphenylmethane (90 mg, $360\mu\text{mol}$, 10 mole equivalents) in NMP. The mixture is heated at 80°C for 16 hours. After being


cooled to room temperature, the reagent solution is decanted and the Lanterns washed with DMF ($3\times 3\text{min}$) and DCM ($3\times 3\text{min}$).



Quaternization

Each D-series Lantern is treated with 0.5 mL of a solution of allyl bromide ($19\mu\text{L}$, $180\mu\text{mol}$, 5 mole

equivalents) in DMF at rt for 16 hours. The reagent solution is decanted and the Lanterns washed with DMF ($3\times 3\text{min}$) and DCM ($3\times 3\text{min}$).




Cleavage (*Hofmann elimination*)

Method A: using DIEA as the base

Each **D-Series Lantern** is treated with 0.5 mL of a solution of DIEA (15 μ L, 72 μ mol, 2 mole equivalents) in DCM at rt for 16 hours. The Lantern is removed (**caution!** retain the reaction solution) and washed with DCM (3 \times 3min). The reaction solution is combined with the DCM washing solutions. The combined solution is concentrated.

The resultant residue is loaded on to a column of


silica gel (50mg, Merck Silica gel 60, 230-400 mesh) and eluted with hexane (3mL), ethyl acetate (5mL) and 20% MeOH/ethyl acetate (5mL) to yield 3.2mg of product. The overall yield is 31%⁵, based on the initial loading of hydroxymethyl Lanterns. Samples are dissolved in 90% CH₃CN/H₂O and in CDCl₃ for LC-MS and ¹H NMR analysis respectively.



Method B: using aminomethyl Lanterns as the base



Each **D-Series Lantern** is treated with 2mL of DCM and 2 PS aminomethyl D-Series Lanterns⁴ (72 μ mol, 2 mole equivalents) at rt for 16 hours with gentle agitation. The Lanterns are removed (**caution!** retain the reaction solution) and washed with DCM (3 \times 3min). The reaction solution is combined with the DCM washing solutions. The

combined solution is concentrated to yield 4.0mg of the final product. The overall yield is 38%⁵, based on the initial loading of hydroxymethyl Lanterns. Samples are dissolved in 90% CH₃CN/H₂O and in CDCl₃ for LC-MS and ¹H NMR analysis respectively.



Reverse phase HPLC trace of the tertiary amine product (Method A and B)


Detection at 214nm



Electrospray MS trace of LC peak at $R_t = 6.27$ min (Method A and B)

Molecular Formula: $C_{20}H_{24}N_2$
 Monoisotopic Mol. Weight: 292.43
 $[M+H]^+$ peak at 293.2amu

400MHz 1H NMR spectrum of the tertiary amine product obtained via Method B ($CDCl_3$)



Notes and References

- 1 Heinonen, P. and Lonnberg, H., *Tetrahedron Lett.*, 1997, **38**, 8569-8572.
- 2 Brown, A.R., Rees, D.C., Rankovic, Z. and Porphy, J.R., *J. Am. Chem. Soc.*, 1997, **119**, 3288-3295.
- 3 Ouyang, X., Armstrong, R.W. and Murphy, M.M., *J. Org. Chem.*, 1998, **63**, 1027-1032.
- 4 PS aminomethyl D-Series Lanterns supplied as TFA salt are neutralised with 5% DIEA in DMF/DCM (v/v, 1:1) (2×10 min), then washed with DMF(3×3 min) and DCM(3×3 min), air dried.
- 5 Yields are comparable to those reported in references 1 and 2.



International
 Tel: + 61 3 9565 1111
 Fax: + 61 3 9565 1199
 mimotopes@mimotopes.com

France
 Tel: + 33 1 5858 0002
 Fax: + 33 1 5858 0006
 europe@mimotopes.com

United Kingdom
 Tel: +44 151 648 3343
 Fax: +44 151 648 3328
 uk@mimotopes.com

USA West
 Tel: + 1 858 558 5800
 Fax: + 1 858 558 5810
 Tel: 800 644 1866
 Fax: 800 655 1866
 uswest@mimotopes.com

USA East
 Tel: + 1 919 873 1123
 Fax: + 1 919 873 1127
 Tel: 800 633 8161
 Fax: 800 424 3970
 useast@mimotopes.com