



Peptide Synthesis Protocol on SynPhase™ PA Lanterns

The **hydrophilic surface** of SynPhase PA (polyamide) Lanterns is based on a proprietary grafted dimethylacrylamide co-polymer, and is efficiently wetted in both polar (including aqueous) and non polar solvents. Two distinct loading levels in the PA range cater for different applications;

- ▶ a low loading version is optimized for peptide synthesis.
- ▶ a high loading version is ideal for small organic molecule synthesis, particularly in polar solvents when polystyrene (PS) is not suitable.

A general synthetic protocol for a complete amino acid coupling cycle using low loading Rink Amide PA Lanterns is described below. The procedure benefits from the convenient handling and washing techniques derived from the modular nature of SynPhase Lanterns, and is readily applied to multiple peptide synthesis as previously described for the SynPhase product range¹. Furthermore, SynPhase PA Lanterns can be readily used in commercial peptide synthesizers². The modular format of the Lanterns means that no weighing is required.

Materials and Methods

Peptides can be synthesized on SynPhase PA RAM D-series Lanterns (loading 8μmol/Lantern, catalog code: SP-PA-D-RAM-008).

Fmoc protected amino acids should be used throughout.

Side chain protection can be afforded by: Trt for Cysteine, Histidine, Glutamine and Asparagine; Boc for Lysine and Tryptophan; tBu for Tyrosine, Threonine, Serine, Aspartic acid and Glutamic acid; Pbf for Arginine.

The general coupling conditions for SynPhase

PA Lanterns utilizes a solution of 80% distilled DMF and 20% AR grade DCM, with reagents at the following ratios and concentrations:

	AA	HOBt	DIC
Reagent ratio	1	1.2	1
Final conc. (mM)	120	144	120

This activated amino acid solution is prepared from equal volumes of an amino acid stock solution and activator solution. For each Lantern a minimum of 500μL of the activated amino acid solution is required.

AA:	Fmoc-amino acid-OH	HOBt:	N-hydroxybenzotriazole (monohydrate)
AR:	analytical reagent	Pbf:	2,2,4,6,7-pentamethylidihydrobenzofuran-5-sulfonyl
Boc:	t-butoxycarbonyl	pip:	piperidine
DCM:	dichloromethane	rt:	room temperature
DIC:	N,N'-diisopropylcarbodiimide	tBu:	t-butyl
DMF:	dimethylformamide	TFA:	trifluoroacetic acid
EDT:	ethanedithiol	TFFH:	tetramethylfluoroformamidinium hexafluorophosphate
Fmoc:	9-fluorenylmethoxycarbonyl	Trt:	trityl
HATU:	N-[(dimethylamino)-1H-1,2,3-triazolo[4,5-b]pyridin-1-ylmethylene]-N-methylmethan-aminium hexafluorophosphate N-oxide		

Reagent Preparation

Amino acid stock solution:

Prepare a 240mM amino AA solution (240mM in 60% distilled DMF and 40% AR grade DCM).

Activator solution:


Prepare a solution containing 288mM HOBt and 240mM of DIC in distilled DMF. An amount equal

to the AA stock solution volume being used for each AA used in the coupling cycle is required. The activator solution should be prepared fresh and used immediately.

Fmoc Group Deprotection solution:

Prepare a solution of 20% piperidine in distilled DMF.

Coupling Scheme



SynPhase D-series Lanterns require typical working volumes of 500 μ L. Actual volumes depend upon the dimensions of the reaction vessel and the number of Lanterns used. The procedure may require slight optimization

depending upon the nature of the peptide(s). For difficult couplings, two hour coupling times and/or the use of sonication may be beneficial. The use of HATU or TFFH is recommended for sterically demanding couplings.

Side Chain Deprotection

Simultaneous side chain deprotection and cleavage can be carried out using 2.5mL per Lantern of a solution of 82.5% TFA/5% thioanisole/5% anisole/5% water/2.5% EDT for 2hr. The TFA solution can be reduced under

vacuum (or a gentle stream of N₂) and the cleaved peptide precipitated in 8mL cold diethyl ether/petroleum ether (bp 40-60°C) (1:2v/v), washed with 4mL cold diethylether/petroleum ether, and then air dried.

References

- 1 Maeji, N.J., Bray, A.M., Valerio, R.M. and Wang, W., Larger Scale Multipin™ Peptide Synthesis, *Peptide Res.*, 1995, **8**, 33-38.
- 2 See SynPhase Application Note SAN002, Automated Fmoc-Peptide Synthesis: Use of SynPhase™ PA Lanterns and PS-Resin in an ACT 496 Ω Multiple Synthesizer.



International
Tel: + 61 3 9565 1111
Fax: + 61 3 9565 1199
mimotopes@mimotopes.com

France
Tel: + 33 1 5858 0002
Fax: + 33 1 5858 0006
europe@mimotopes.com

United Kingdom
Tel: + 44 151 648 3343
Fax: + 44 151 648 3328
uk@mimotopes.com

USA West
Tel: + 1 858 558 5800
Fax: + 1 858 558 5810
Tel: 800 644 1866
Fax: 800 655 1866
uswest@mimotopes.com

USA East
Tel: + 1 919 873 1123
Fax: + 1 919 873 1127
Tel: 800 633 8161
Fax: 800 424 3970
useast@mimotopes.com