




Solid Phase Synthesis of Benzimidazoles: A Comparative Study on SynPhase™ A-, D- and L-series Lanterns

SynPhase A-series Lanterns are the latest addition to the SynPhase product range and offer a higher loading than L- and D-series Lanterns. All three series offer the same performance in terms of kinetics, yield and purity, and chemistry transfers between the different Lantern types without the need for re-optimization.

This technical note demonstrates the ease of chemistry transfer using a one-pot reduction-cyclization synthesis of benzimidazoles, previously developed by Mimotopes on SynPhase PS D-series Lanterns¹.

Materials and Methods



DCM: dichloromethane
DIC: *N,N'*-diisopropylcarbodiimide
DMF: dimethylformamide
Fmoc: 9-fluorenylmethoxycarbonyl
HOBT: *N*-hydroxybenzotriazole (monohydrate)

pip: piperidine
DIEA: diisopropylethylamine
TFA: trifluoroacetic acid
NMP: 1-methyl-2-pyrrolidinone

Reagent Preparation

SynPhase PS Fmoc Rink amide L-, D- and A-series Lanterns (15 each, loading 15, 35, 75 μ mol respectively) in a reaction flask were treated with 20% piperidine (30mL) at room temperature for 45min. The reagent solution was decanted and Lanterns washed with DMF (3 \times 3min), DCM (3 \times 3min) and air dried.

The Fmoc deprotected Lanterns were treated with 30mL of a solution of 4-fluoro-3-nitrobenzoic acid (0.15M), DIC (0.15M) and HOBt (0.075M) in 20% DMF/DCM at room temperature for 16h. The reagent solution was decanted and Lanterns washed with DMF (3 \times 3min), DCM (3 \times 3min) and air dried to give the Lantern-bound **2**.

Lanterns **2** were treated with 30mL of a solution of n-propylamine (0.5M) in 5% DIEA/NMP at 60°C

for 16h. The solution was decanted and Lanterns washed with DMF (3 \times 3min), DCM (3 \times 3min) and air dried to give the Lantern-bound **3**.


Four Lanterns each of L-, D- and A-series **3** were treated with (1.0, 2.0 & 4.0mL respectively) of a solution of tin(II) chloride dihydrate (0.75M) and 4-(trifluoromethyl)benzaldehyde (0.14M, 2 eq) in DMF at 60°C for 16h. The reagent solution was decanted and Lanterns washed with DMF (3 \times 3min), DCM (3 \times 3min) and air dried to give the Lantern-bound **4** and **5**. Each group of 4 L-, D- and A-series Lanterns was cleaved with 20% TFA/DCM (2.0, 3.0 & 4.0 mL respectively) at room temperature for 1h. The cleavage solution was evaporated to dryness. The residue was dissolved in 10% H₂O/CH₃CN for HPLC and LC-MS analysis.

Results and Conclusion

All benzimidazole products **6** obtained from three different types of Lanterns gave nearly identical HPLC chromatograms and purities. The purity of the products were 76 (19)%, 78 (17)% and 74 (19)% corresponding to L-, D- and A-series Lanterns respectively. The figures in parentheses indicate the percentage of the by-product benzimidazolium **7**. Both the benzimidazole **6** and benzimidazolium **7** are confirmed by LC-MS. The yield was >95% in all cases.

The one-pot benzimidazole synthesis on SynPhase Lanterns has demonstrated that the chemistry developed on a particular series of Lanterns (D-series in this case) can be freely

transferred between L-, D- and A-series Lanterns without compromise of results.



Reverse phase HPLC trace of a crude benzimidazole product synthesized on A-, D- and L-series Lanterns under identical conditions.

Detection at 214nm

References

- 1 Wu, Z., Rea, P., Wickham, G., Tetrahedron Letters, 2000, **41**, 9871-9874.



International
Tel: + 61 3 9565 1111
Fax: + 61 3 9565 1199
mimotopes@mimotopes.com

France
Tel: + 33 1 5858 0002
Fax: + 33 1 5858 0006
europe@mimotopes.com

United Kingdom
Tel: +44 151 648 3343
Fax: +44 151 648 3328
uk@mimotopes.com

USA West
Tel: + 1 858 558 5800
Fax: + 1 858 558 5810
Tel: 800 644 1866
Fax: 800 655 1866
uswest@mimotopes.com

USA East
Tel: + 1 919 873 1123
Fax: + 1 919 873 1127
Tel: 800 633 8161
Fax: 800 424 3970
useast@mimotopes.com